


Product Datasheet

53BP1 (TP53BP1) antibody (orb1316148)

Description	Mouse monoclonal antibody to 53BP1 (TP53BP1)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-300 of human TP53BP1(NP_005648) produced in E.coli.
Target	TP53BP1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	213.4 kDa
Uniprot ID	Q12888
NCBI	NM_005657
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...