


Product Datasheet

Caspase 4 (CASP4) antibody (orb1316819)

| | |
|----------------------------|---|
| Description | Mouse monoclonal antibody to Caspase 4 (CASP4) |
| Species/Host | Mouse |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Immunogen | Human recombinant protein fragment corresponding to amino acids 81-270 of human CASP4 (NP_150649) produced in E.coli. |
| Target | CASP4 |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| MW | 36.5 kDa |
| Uniprot ID | P49662 |
| NCBI | NM_033306 |
| Dilution Range | WB 1:2000 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of HEK293T cells w...