


Product Datasheet

hnRNP L (HNRNPL) antibody (orb1316895)

| | |
|----------------------------|---|
| Description | Mouse monoclonal antibody to hnRNP L (HNRNPL) |
| Species/Host | Mouse |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | IHC, WB |
| Immunogen | Human recombinant protein fragment corresponding to amino acids 270-502 of human HNRNPL (NP_001524) produced in E.coli. |
| Target | HNRNPL |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| MW | 64 kDa |
| Uniprot ID | P14866 |
| NCBI | NM_001533 |
| Dilution Range | WB 1:2000, IHC 1:150 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of HEK293T cells w...