

Product Datasheet

CALD1 antibody pairs (orb131724)

Description	Rabbit polyclonal to CALD1. CALD1 antibodies represent a pair of fully characterized antibodies that...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal or C-terminal region of human CALD1
Preservatives	0.01% NaN ₃
Concentration	100 µg/ml BSA
Storage	-20°C
Note	For research use only
Application notes	<p>WB: The apparent protein size on WB may be different from the calculated M.W. due to modifications.</p> <p>Experiment Notes: Rabbit polyclonal antibodies were produced by immunizing animals with GST-fusion proteins containing either the N-terminal or the C-terminal region of human CALD1.</p>
Clonality	Polyclonal
Purity	Affinity purification
MW	93
Uniprot ID	Q05682
Dilution Range	WB:1:1000-1:3000
Expiration Date	12 months from date of receipt.