

## Product Datasheet

### CALD1 antibody (orb131726)

<b>Description</b>	Rabbit polyclonal to CALD1. Caldesmon (CALD1) belongs to the caldesmon family and is an actin-and my...
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	C-terminal region of human CALD1
<b>Preservatives</b>	0.01% NaN <sub>3</sub>
<b>Concentration</b>	100 µg/ml BSA
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Application notes</b>	<p><b>WB:</b> The apparent protein size on WB may be different from the calculated M.W. due to modifications.</p> <p><b>Experiment Notes:</b> Rabbit polyclonal antibodies were produced by immunizing animals with a GST-fusion protein containing the C-terminal region of human CALD1.</p>
<b>Clonality</b>	Polyclonal
<b>Purity</b>	Affinity purification
<b>MW</b>	93
<b>Uniprot ID</b>	<a href="#">Q05682</a>
<b>Dilution Range</b>	WB:1:1000-1:3000
<b>Expiration Date</b>	12 months from date of receipt.