


## Product Datasheet

### CD39 (ENTPD1) antibody (orb1317289)

<b>Description</b>	Mouse monoclonal antibody to CD39 (ENTPD1)
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC, WB
<b>Immunogen</b>	Human recombinant protein fragment corresponding to amino acids 38-307 of human ENTPD1 (NP_001767) produced in E.coli.
<b>Target</b>	ENTPD1
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<a href="#">P49961</a>
<b>NCBI</b>	<a href="#">NM_001776</a>
<b>Dilution Range</b>	WB 1:2000, IHC 1:150
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of Human recombinant protein fragment corresponding to amino acids 38-307 of human ENTPD1 (NP\_001767) produced in E.coli.