


Product Datasheet

SP17 (SPA17) antibody (orb1317407)

Description	Mouse monoclonal antibody to SP17 (SPA17)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human SPA17 (NP_059121) produced in E.coli.
Target	SPA17
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	17.2 kDa
Uniprot ID	Q15506
NCBI	NM_017425
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...