

## Product Datasheet

### Calretinin (CALB2) antibody (orb1317428)

<b>Description</b>	Mouse monoclonal antibody to Calretinin (CALB2)
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	Full length human recombinant protein of human CALB2 (NP_001731) produced in HEK293T cell.
<b>Target</b>	CALB2
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>MW</b>	31.4 kDa
<b>Uniprot ID</b>	<a href="#">P22676</a>
<b>NCBI</b>	<a href="#">NM_001740</a>
<b>Dilution Range</b>	WB 1:500
<b>Expiration Date</b>	12 months from date of receipt.

