

Product Datasheet

UBE3A antibody (orb131783)

Description Rabbit polyclonal to UBE3A. Ubiquitin-protein ligase E3A (UBE3A) is an E3 ubiquitin-protein ligase w...

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications WB

Immunogen N-terminal region of human UBE3A

Preservatives 0.01% NaN₃

Concentration 100 µg/ml BSA

Storage -20°C

Note For research use only

Application notes **WB:** The apparent protein size on WB may be different from the calculated M.W. due to modifications.

Experiment Notes: Rabbit polyclonal antibodies were produced by immunizing animals with a GST-fusion protein containing the N-terminal region of human UBE3A. Antibodies were purified by affinity purification using immunogen.

Clonality Polyclonal

Purity Affinity purification

MW 101

Uniprot ID **Q05086**

Dilution Range WB:1:1000-1:3000

Expiration Date 12 months from date of receipt.