


Product Datasheet

Galectin 2 (LGALS2) antibody (orb1317924)

Description	Mouse monoclonal antibody to Galectin 2 (LGALS2)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human LGALS2 (NP_006489) produced in E.coli.
Target	LGALS2
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	14.5 kDa
Uniprot ID	P05162
NCBI	NM_006498
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...