

Product Datasheet

SETDB2 antibody (orb131820)

Description	Rabbit polyclonal to SETDB2. Histone-lysine N-methyltransferase SETDB2 belongs to the histone-lysine...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human SETDB2
Preservatives	0.01% NaN ₃
Concentration	100 µg/ml BSA
Storage	-20°C
Note	For research use only
Application notes	<p>WB: The apparent protein size on WB may be different from the calculated M.W. due to modifications.</p> <p>Experiment Notes: Rabbit polyclonal antibodies were produced by immunizing animals with a GST-fusion protein containing the N-terminal region of human SETDB2.</p>
Clonality	Polyclonal
Purity	Affinity purification
MW	82
Uniprot ID	Q96T68
Dilution Range	WB:1:1000-1:3000
Expiration Date	12 months from date of receipt.