


Product Datasheet

VILIP1 (VSNL1) antibody (orb1318475)

Description	Mouse monoclonal antibody to VILIP1 (VSNL1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 2-191 of human VSNL1 (NP_003376) produced in E.coli.
Target	VSNL1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	22 kDa
Uniprot ID	P62760
NCBI	NM_003385
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...