


## Product Datasheet

### VILIP1 (VSNL1) antibody (orb1318478)

<b>Description</b>	Mouse monoclonal antibody to VILIP1 (VSNL1)
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC, WB
<b>Immunogen</b>	Human recombinant protein fragment corresponding to amino acids 2-191 of human VSNL1 (NP_003376) produced in E.coli.
<b>Target</b>	VSNL1
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>MW</b>	22 kDa
<b>Uniprot ID</b>	<a href="#">P62760</a>
<b>NCBI</b>	<a href="#">NM_003385</a>
<b>Dilution Range</b>	WB 1:2000, IHC 1:150
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of HEK293T cells w...