


## Product Datasheet

# Thymidylate Synthase (TYMS) antibody (orb1318608)

<b>Description</b>	Mouse monoclonal antibody to Thymidylate Synthase (TYMS)
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	Full length human recombinant protein of human TYMS (NP_001062) produced in E.coli.
<b>Target</b>	TYMS
<b>Note</b>	For research use only
<b>Isotype</b>	IgG2a
<b>Clonality</b>	Monoclonal
<b>MW</b>	35.5 kDa
<b>Uniprot ID</b>	<a href="#">P04818</a>
<b>NCBI</b>	<a href="#">NM_001071</a>
<b>Dilution Range</b>	WB 1:2000
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of HEK293T cells w...