

Product Datasheet

NPAS3 antibody (orb131865)

Description	Rabbit polyclonal to NPAS3. Neuronal PAS domain-containing protein 3 (NPAS3) is a member of the basi...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human NPAS3
Preservatives	0.01% NaN ₃
Concentration	100 µg/ml BSA
Storage	-20°C
Note	For research use only
Application notes	<p>WB: The apparent protein size on WB may be different from the calculated M.W. due to modifications.</p> <p>Experiment Notes: Rabbit polyclonal antibodies were produced by immunizing animals with a GST-fusion protein containing the C-terminal region of human NPAS3.</p>
Clonality	Polyclonal
Purity	Affinity purification
MW	101
Uniprot ID	Q8IXF0
Dilution Range	WB:1:1000-1:3000
Expiration Date	12 months from date of receipt.