

Product Datasheet

PPRC1 antibody (orb131866)

Description	Rabbit polyclonal to PPRC1. Peroxisome proliferator-activated receptor gamma coactivator-related pro...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal region of human PPRC1
Preservatives	0.01% NaN ₃
Concentration	100 µg/ml BSA
Storage	-20°C
Note	For research use only
Application notes	<p>WB: The apparent protein size on WB may be different from the calculated M.W. due to modifications.</p> <p>Experiment Notes: Rabbit polyclonal antibodies were produced by immunizing animals with a GST-fusion protein containing the N-terminal region of human PPRC1.</p>
Clonality	Polyclonal
Purity	Affinity purification
MW	178
Uniprot ID	Q5VV67
Dilution Range	WB:1:1000-1:3000
Expiration Date	12 months from date of receipt.