


## Product Datasheet

### CD23 (FCER2) antibody (orb1318686)

<b>Description</b>	Mouse monoclonal antibody to CD23 (FCER2)
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	FC, WB
<b>Immunogen</b>	Human recombinant protein fragment corresponding to amino acids 48-321 of human FCER2 (NP_001993) produced in SF9 cell.
<b>Target</b>	FCER2
<b>Note</b>	For research use only
<b>Isotype</b>	IgG2b
<b>Clonality</b>	Monoclonal
<b>MW</b>	36.3 kDa
<b>Uniprot ID</b>	<a href="#">P06734</a>
<b>NCBI</b>	<a href="#">NM_002002</a>
<b>Dilution Range</b>	WB 1:2000
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of HEK293T cells w...