

Product Datasheet

TJP3 antibody (orb131875)

| | |
|----------------------------|--|
| Description | Rabbit polyclonal to TJP3. Tight junction protein 3 (TJP3) is a member of the membrane-associated gu... |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Immunogen | C-terminal region of human TJP3 |
| Preservatives | 0.01% NaN ₃ |
| Concentration | 100 µg/ml BSA |
| Storage | -20°C |
| Note | For research use only |
| Application notes | <p>WB: The apparent protein size on WB may be different from the calculated M.W. due to modifications.</p> <p>Experiment Notes: Rabbit polyclonal antibodies were produced by immunizing animals with a GST-fusion protein containing the C-terminal region of human TJP3.</p> |
| Clonality | Polyclonal |
| Purity | Affinity purification |
| MW | 103 |
| Uniprot ID | O95049 |
| Dilution Range | WB:1:1000-1:3000 |
| Expiration Date | 12 months from date of receipt. |