

Product Datasheet

EFNA5 antibody (orb131880)

Description	Rabbit polyclonal to EFNA5. Ephrin-A5, a member of the ephrin gene family, prevents axon bundling in...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human EFNA5
Preservatives	0.01% NaN ₃
Concentration	100 µg/ml BSA
Storage	-20°C
Note	For research use only
Application notes	<p>WB: The apparent protein size on WB may be different from the calculated M.W. due to modifications.</p> <p>Experiment Notes: Rabbit polyclonal antibodies were produced by immunizing animals with a GST-fusion protein containing the C-terminal region of human EFNA5.</p>
Clonality	Polyclonal
Purity	Affinity purification
MW	26
Uniprot ID	P52803
Dilution Range	WB:1:1000-1:3000
Expiration Date	12 months from date of receipt.