


## Product Datasheet

# Isocitrate dehydrogenase (IDH1) antibody (orb1319267)

<b>Description</b>	Mouse monoclonal antibody to Isocitrate dehydrogenase (IDH1)
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Canine, Human, Monkey, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC, WB
<b>Immunogen</b>	Synthetic peptide around the R132 region of the human IDH conjugated to KLH
<b>Target</b>	IDH1
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG2b
<b>Clonality</b>	Monoclonal
<b>MW</b>	46.5 kDa
<b>Uniprot ID</b>	<a href="#">O75874</a>
<b>NCBI</b>	<a href="#">NM_005896</a>
<b>Dilution Range</b>	WB 1:2000, IHC 1:150
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of Equivalent amou...