


Product Datasheet

Isocitrate dehydrogenase (IDH1) antibody (orb1319280)

Description	Mouse monoclonal antibody to Isocitrate dehydrogenase (IDH1)
Species/Host	Mouse
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Synthetic peptide around the R132H mutation region of the human IDH conjugated to KLH
Target	IDH1
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	46.5 kDa
Uniprot ID	O75874
NCBI	NM_005896
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...