


Product Datasheet

B Raf (BRAF) antibody (orb1319590)

Description	Rabbit monoclonal antibody to B Raf (BRAF)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic peptide (the amino acid sequence is considered to be commercially sensitive) within Human BRAF V600E (NP_004324). The exact sequence is proprietary.
Target	BRAF
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P15056
NCBI	NM_004333
Dilution Range	WB1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of overexpressed I...