


Product Datasheet

His tag antibody (orb1319747)

Description	Rabbit monoclonal antibody to His tag
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic peptide conjugated to KLH
Target	His
Preservatives	Tissue culture buffer containing 1% rAlbumin and 0.02% Sodium azide
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Dilution Range	WB: 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from S...