


Product Datasheet

Neurofibromin (NF1) antibody (orb1319813)

Description	Rabbit polyclonal antibody to Neurofibromin (NF1)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	DNA immunization. This antibody is specific for the N Terminus Region of the target protein.
Target	NF1
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P21359
NCBI	NM_001042492
Dilution Range	ELISA: 1:100-1:2000
Expiration Date	12 months from date of receipt.



ELISA using Neurofibromin (NF1) antibody