


Product Datasheet

Moesin (MSN) antibody (orb1319842)

Description	Rabbit polyclonal antibody to Moesin (MSN)
Species/Host	Rabbit
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	DNA immunization. This antibody is specific for the C Terminus Region of the target protein.
Target	MSN
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P26038
NCBI	NM_002444
Dilution Range	WB 1:5000~20000, IF 1:500,ELISA 1:100-1:2000
Expiration Date	12 months from date of receipt.



ELISA using Moesin (MSN) antibody