


Product Datasheet

PTP1B (PTPN1) antibody (orb1319850)

Description	Rabbit polyclonal antibody to PTP1B (PTPN1)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	DNA immunization. In vivo generated recombinant protein fragment
Target	PTPN1
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P18031
NCBI	NM_002827
Dilution Range	WB 1:20000, IHC 1:150, ELISA 1:100-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...