


Product Datasheet

Liprin alpha 1 (PPFIA1) antibody (orb1320024)

Description	Rabbit polyclonal antibody to Liprin alpha 1 (PPFIA1)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	DNA immunization. This antibody is specific for the Middle Region of the target protein.
Target	PPFIA1
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q13136
NCBI	NM_177423
Dilution Range	WB: 1:5000-1:20000; ELISA: 1:100-1:2000; IHC: 1:10-1:2000; IHC-P 1:10-1:2000
Expiration Date	12 months from date of receipt.



Western
Blot
analysis of
Liprin alpha
1 ...