


## Product Datasheet

### Nucleophosmin (NPM1) antibody (orb1320040)

<b>Description</b>	Rabbit polyclonal antibody to Nucleophosmin (NPM1)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	DNA immunization. This antibody is specific for the C Terminus Region of the target protein.
<b>Target</b>	NPM1
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<a href="#">P06748</a>
<b>NCBI</b>	<a href="#">NM_199185</a>
<b>Dilution Range</b>	WB: 1:5000-1:20000; ELISA: 1:100-1:2000; IHC: 1:10-1:2000; IHC-P 1:250-1:2000
<b>Expiration Date</b>	12 months from date of receipt.



ELISA using Nucleophosmin (NPM1) antibody...