


## Product Datasheet

### Prion protein PrP (PRNP) antibody (orb1320061)

<b>Description</b>	Rabbit polyclonal antibody to Prion protein PrP (PRNP)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	DNA immunization. This antibody is specific for the N Terminus Region of the target protein.
<b>Target</b>	PRNP
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<a href="#">F7VJQ1</a>
<b>NCBI</b>	<a href="#">NM_183079</a>
<b>Dilution Range</b>	ELISA: 1:100-1:2000
<b>Expiration Date</b>	12 months from date of receipt.



ELISA using Prion protein PrP (PRNP) ant...