


Product Datasheet

MFI2 (MELTF) antibody (orb1320116)

Description	Rabbit polyclonal antibody to MFI2 (MELTF)
Species/Host	Rabbit
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	DNA immunization. This antibody is specific for the N Terminus Region of the target protein.
Target	MELTF
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P08582
NCBI	NM_005929
Dilution Range	WB 1:5000~20000, ELISA 1:100-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...