


Product Datasheet

Munc 13 4 (UNC13D) antibody (orb1320279)

Description	Mouse monoclonal antibody to Munc 13 4 (UNC13D)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human UNC13D(NP_954712) produced in HEK293T cell.
Target	UNC13D
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	123.1 kDa
Uniprot ID	Q70J99
NCBI	NM_199242
Dilution Range	WB 1:200~1000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...