


Product Datasheet

SCF (KITLG) antibody (orb1320310)

Description	Mouse monoclonal antibody to SCF (KITLG)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human KITLG(NP_000890) produced in HEK293T cell.
Target	KITLG
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	27.9 kDa
Uniprot ID	P21583
NCBI	NM_000899
Dilution Range	WB 1:4000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...