


## Product Datasheet

### SIVA (SIVA1) antibody (orb1320525)

<b>Description</b>	Mouse monoclonal antibody to SIVA (SIVA1)
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IF, IHC, WB
<b>Immunogen</b>	Full length human recombinant protein of human SIVA1(NP_006418) produced in HEK293T cell.
<b>Target</b>	SIVA1
<b>Note</b>	For research use only
<b>Isotype</b>	IgG2a
<b>Clonality</b>	Monoclonal
<b>MW</b>	18.5 kDa
<b>Uniprot ID</b>	<a href="#">O15304</a>
<b>NCBI</b>	<a href="#">NM_006427</a>
<b>Dilution Range</b>	WB 1:2000, IF: 1:100
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of HEK293T cells w...