


Product Datasheet

SIVA (SIVA1) antibody (orb1320526)

Description	Mouse monoclonal antibody to SIVA (SIVA1)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human SIVA1(NP_006418) produced in HEK293T cell.
Target	SIVA1
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	18.5 kDa
Uniprot ID	O15304
NCBI	NM_006427
Dilution Range	WB 1:1000, IF: 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...