


Product Datasheet

alpha Synuclein (SNCA) antibody (orb1320655)

Description	Mouse monoclonal antibody to alpha Synuclein (SNCA)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human SNCA(NP_000336) produced in HEK293T cell.
Target	SNCA
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	14.3 kDa
Uniprot ID	P37840
NCBI	NM_000345
Dilution Range	WB 1:4000, IHC: 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...