


Product Datasheet

Lysozyme (LYZ) antibody (orb1320766)

Description	Mouse monoclonal antibody to Lysozyme (LYZ)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full length human recombinant protein of human LYZ(NP_000230) produced in HEK293T cell.
Target	LYZ
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	14.7 kDa
Uniprot ID	P61626
NCBI	NM_000239
Dilution Range	WB 1:4000, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...