


Product Datasheet

IL36 gamma (IL36G) antibody (orb1320915)

Description	Mouse monoclonal antibody to IL36 gamma (IL36G)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human IL1F9(NP_062564) produced in HEK293T cell.
Target	IL36G
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	18.5 kDa
Uniprot ID	Q9NZH8
NCBI	NM_019618
Dilution Range	WB 1:4000, IHC 1:150, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...