


Product Datasheet

Hsc70 (HSPA8) antibody (orb1320936)

Description	Mouse monoclonal antibody to Hsc70 (HSPA8)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human HSPA8(NP_006588) produced in HEK293T cell.
Target	HSPA8
Note	For research use only
Isotype	IgG3
Clonality	Monoclonal
MW	70.7 kDa
Uniprot ID	P11142
NCBI	NM_006597
Dilution Range	WB 1:1000, IHC 1:150, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...