


Product Datasheet

MAPK4 antibody (orb1321005)

Description	Mouse monoclonal antibody to MAPK4
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human MAPK4(NP_002738) produced in HEK293T cell.
Target	MAPK4
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	65.7 kDa
Uniprot ID	P31152
NCBI	NM_002747
Dilution Range	WB 1:4000, IHC 1:150, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...