


## Product Datasheet

### ERp57 (PDIA3) antibody (orb1321679)

<b>Description</b>	Mouse monoclonal antibody to ERp57 (PDIA3)
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC, WB
<b>Immunogen</b>	Human recombinant protein fragment corresponding to amino acids 140-505 of human PDIA3(NP_005304) produced in E.coli.
<b>Target</b>	PDIA3
<b>Note</b>	For research use only
<b>Isotype</b>	IgG2a
<b>Clonality</b>	Monoclonal
<b>MW</b>	54.2 kDa
<b>Uniprot ID</b>	<a href="#">P30101</a>
<b>NCBI</b>	<a href="#">NM_005313</a>
<b>Dilution Range</b>	WB 1:1000, IHC 1:150
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of HEK293T cells w...