


Product Datasheet

ERp57 (PDIA3) antibody (orb1321686)

Description	Mouse monoclonal antibody to ERp57 (PDIA3)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids PDIA3(NP_005304) produced in E.coli.
Target	PDIA3
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	54.2 kDa
Uniprot ID	P30101
NCBI	NM_005313
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...