


Product Datasheet

Annexin VII (ANXA7) antibody (orb1321734)

Description	Mouse monoclonal antibody to Annexin VII (ANXA7)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 105-411 of human ANXA7(NP_004025) produced in E.coli.
Target	ANXA7
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	52.6 kDa
Uniprot ID	P20073
NCBI	NM_004034
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...