


## Product Datasheet

### JAB1 (COPS5) antibody (orb1322163)

<b>Description</b>	Mouse monoclonal antibody to JAB1 (COPS5)
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IF, IHC, WB
<b>Immunogen</b>	Full length human recombinant protein of human COPS5(NP_006828) produced in HEK293T cell.
<b>Target</b>	COPS5
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>MW</b>	37.4 kDa
<b>Uniprot ID</b>	<a href="#">Q92905</a>
<b>NCBI</b>	<a href="#">NM_006837</a>
<b>Dilution Range</b>	WB 1:2000, IHC 1:150, IF 1:100
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of HEK293T cells w...