


## Product Datasheet

# Cytochrome P450 2B6 (CYP2B6) antibody (orb1322178)

<b>Description</b>	Mouse monoclonal antibody to Cytochrome P450 2B6 (CYP2B6)
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Canine, Human, Monkey, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	FC, IF, IHC, WB
<b>Immunogen</b>	Full length human recombinant protein of human CYP2B6(NP_000758) produced in HEK293T cell.
<b>Target</b>	CYP2B6
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG2a
<b>Clonality</b>	Monoclonal
<b>MW</b>	56.1 kDa
<b>Uniprot ID</b>	<a href="#">P20813</a>
<b>NCBI</b>	<a href="#">NM_000767</a>
<b>Dilution Range</b>	WB 1:500~2000, IHC 1:150, IF 1:100, FLOW 1:100
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of HEK293T cells w...