


Product Datasheet

C5orf19 (REEP2) antibody (orb1322182)

Description	Mouse monoclonal antibody to C5orf19 (REEP2)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human REEP2(NP_057690) produced in HEK293T cell.
Target	REEP2
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	28.1 kDa
Uniprot ID	Q9BRK0
NCBI	NM_016606
Dilution Range	WB 1:500~2000, IHC 1:150, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...