


## Product Datasheet

# Coatomer subunit delta (ARCN1) antibody (orb1322187)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Mouse monoclonal antibody to Coatomer subunit delta (ARCN1)   |
| <b>Species/Host</b>        | Mouse   |
| <b>Reactivity</b>          | Canine, Human, Monkey, Mouse, Rat   |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Tested Applications</b> | FC, IF, IHC, WB   |
| <b>Immunogen</b>           | Full length human recombinant protein of human ARCN1(NP_001646) produced in HEK293T cell.   |
| <b>Target</b>              | ARCN1   |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                | For research use only   |
| <b>Isotype</b>             | IgG1  |
| <b>Clonality</b>           | Monoclonal  |
| <b>MW</b>                  | 57 kDa  |
| <b>Uniprot ID</b>          | <a href="#">P48444</a>  |
| <b>NCBI</b>                | <a href="#">NM_001655</a>   |
| <b>Dilution Range</b>      | WB 1:500~2000, IHC 1:50, IF 1:100, FLOW 1:100   |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |



Western blot analysis of HEK293T cells w...