


Product Datasheet

IGJ (JCHAIN) antibody (orb1322293)

Description	Mouse monoclonal antibody to IGJ (JCHAIN)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Full length human recombinant protein of human IGJ(NP_653247) produced in HEK293T cell.
Target	JCHAIN
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	17.9 kDa
Uniprot ID	P01591
NCBI	NM_144646
Dilution Range	WB 1:2000, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...