


## Product Datasheet

### ERp72 (PDIA4) antibody (orb1322331)

<b>Description</b>	Mouse monoclonal antibody to ERp72 (PDIA4)
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Canine, Human, Monkey, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC, WB
<b>Immunogen</b>	Full length human recombinant protei of human PDIA4(NP_004902) produced in HEK293T cell.
<b>Target</b>	PDIA4
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>MW</b>	72.8 kDa
<b>Uniprot ID</b>	<a href="#">P13667</a>
<b>NCBI</b>	<a href="#">NM_004911</a>
<b>Dilution Range</b>	WB 1:500~2000, IHC 1:150
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of HEK293T cells w...