


Product Datasheet

ERp72 (PDIA4) antibody (orb1322518)

Description	Mouse monoclonal antibody to ERp72 (PDIA4)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human PDIA4(NP_004902) produced in HEK293T cell.
Target	PDIA4
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	72.8 kDa
Uniprot ID	P13667
NCBI	NM_004911
Dilution Range	WB 1:2000, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...