


Product Datasheet

TRK fused gene (TFG) antibody (orb1322551)

Description	Mouse monoclonal antibody to TRK fused gene (TFG)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human TFG(NP_001007566) produced in HEK293T cell.
Target	TFG
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	43.3 kDa
Uniprot ID	Q92734
NCBI	NM_001007565
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...